

EXPERIMENT NO.1

NEWTON'S RINGS EXPERIMENT

OBJECT:

To determine the wavelength of Monochromatic (sodium) light by Newton's Ring Method.

APPARATUS:

A plano-convex lens of large radius of curvature, optical arrangement for Newton's Ring, a plane glass plate, a travelling microscope and sodium lamp.

FORMULA USED:

The wavelength (λ) of monochromatic light source is given by

$$\lambda = \frac{D_{n+p}^2 - D_n^2}{4pR}$$

D_{n+p} = Diameter of the $(n+p)^{th}$ ring.

D_n = Diameter of the $(n)^{th}$ ring.

p = An integer number

R = Radius of curvature of the curved surface of the Plano-convex lens.

THEORY:

The Newton's rings are formed by interference between the rays of light reflected from the top and bottom surfaces of the air film formed between the bottom of the Plano-convex lens and the glass plate. The effective path difference between these interfering rays is given as

$$d = 2\mu t \cos r + \frac{\lambda}{2} \quad (1)$$

Where μ is the refractive index of the air film (for air $\mu = 1$); t is thickness of the air film, r is the angle of reflectance and λ is the wavelength of the light.

In the interference pattern, the diameter of the n^{th} bright ring is given by

$$D_n^2 = 2(2n-1)\lambda R \quad (2)$$

Where, R is the radius of curvature of the Plano-convex lens. If D_{n+p} be the diameter of $(n+p)^{th}$ bright ring, then




Figure.1 Experimental Arrangement of Newton's Rings Experiment

$$D_{n+p}^2 = 2[2(n+p)-1]\lambda R$$

Subtracting equation (1) from (2), we get

$$D_{n+p}^2 - D_n^2 = 4p\lambda R$$

or

$$\lambda = \frac{D_{n+p}^2 - D_n^2}{4pR} \text{ A}^0 \quad (3)$$

Equation (3) is used to determine the wavelength of monochromatic light.




Figure.2 Movement of Cross wire




Figure.3. Experimental setup of Newton' Rings Experiment in the Lab

PROCEDURE:

1. Make the experimental arrangement as shown in figure (1) and allow the beam of light to pass through the convex lens L_1 .

[illegible]

CALCULATIONS:

The wavelength of monochromatic light (Sodium lamp) is determined as

$$\lambda = \frac{D_{n+p}^2 - D_n^2}{4pR} = \dots\dots\dots \text{mm} = \dots\dots\dots A^\circ$$

RESULT:

(A) Standard result:

Wavelength of sodium light = $5893 A^\circ$.

(B) Experimental result = $\dots\dots\dots A^\circ$.

(C) % Error = $\dots\dots\dots$.

PRECAUTIONS:

1. The glass plate and lens should be cleaned.
2. The microscope eyepiece must be properly focused on its cross-wire.
3. Microscope should be given motion only along one-direction to avoid back-lash error.
4. The radius of curvature of the surface of the lens in contact with the glass plate should be measured accurately.
5. A lens of large radius of curvature should be used.

REFERENCES:

1. Optics by Brij Lal and Subramaniam.
2. Optics by A.K.Ghatak

VIVA-VOCE

1. What do you understand by the interference of light?
2. What are essential conditions for obtaining interference of light?
3. What do you understand by coherent sources?
4. Is it possible to observe interference pattern by having two independent sources such as two candles?
5. Why should be two sources be monochromatic?
6. Why are the Newton's rings circular?
7. Why is central ring dark?
8. Where are these rings formed?
9. What will happen when a little water is introduced in between the Plano - convex lens and the plate?
10. How does the diameter of rings change on the introduction of liquid?